

WE COULD MAKE THEM THIS WAY...BUT WE DON'T!

Bovie

BOVIE MEDICAL CORPORATION • 5115 Ulmerton Road, Clearwater, FL 33760 • USA • Ph 1-800-537-2790 • sales@boviemed.com • www.boviemed.com



We could make them this way...but we don't! We make them better!

1. We could use heat shrink instead of a double molding process...But we don't.
2. We could leave off the molded ribs for easy insertion and removal from the pencil...But we don't.
3. We could leave out the patented protective tray for the Tungsten loop...But we don't.
4. We could make them in China...But we don't. We make them in the U.S.A.
5. We could skip the 100% inspection, which includes running 3000 volts through the electrode for 60 seconds (double the AAMI/ANSI standard)...But we don't.
6. We could charge a lot more for our loop electrodes...But we don't.



There are many steps we could leave out of our process for making a great loop electrode...but we don't!

LIST OF ELECTRODE STYLES AND CONVERSION CHART...

BOVIE / AARON ELECTRODES		COOPER	PREMIER	WALLACH	CONMED	VALLEY	ACMI/CIRCON	UTAH	MILTEX		
Ball Electrodes	ES06	Extended 3mm Ball, 5/box	B0312	9006225	909038	60-5181-102	E1563	6156901	DBL-311	31-5000	
	ES07	Extended 5mm Ball, 5/box	B0512	9006230	909003	60-5181-103	E1564	6157901	DBL-511	31-5002	
	ES50	1/8" Ball, 25/box									
LLETZ Loop Electrodes	ES08	5mm x 5mm Loop, 5/box						6133901	DLP-T11	31-3056	
	ES09	10mm x 10mm Loop, 5/box	R1010	9006200	909007		E1559	6134901	DLP-S11	31-5072	
	ES10	15mm x 8mm Loop, 5/box	R1507 *	9006205	909011			6135901			
	ES11	15mm x 10mm Loop, 5/box				60-5180-102		6136901	DLP-M11	31-5074	
	ES12	20mm x 8mm Loop, 5/box	R2008	9006215	909013			6137901		31-5090	
	ES13	20mm x 15mm Loop, 5/box	R2008	9006215	909013			6137901		31-5096	
	ES31	20mm x 20mm Loop, 5/box	R2020							31-5098	
	ES42	20mm x 12mm Loop, 5/box		9006220 *				E1561	DLP-W11	31-5094	
	ES43	15mm x 15mm Loop, 5/box	R1515	9006210							
	ES44	15mm x 5mm Loop, 5/box	R1505		909134						
	ES45	13mm x 13mm Loop, 5/box									
	ES46	10mm x 5mm Loop, 5/box									
	ES47	25mm x 10mm Loop, 5/box	R2510								
	ES48	10mm x 3mm Loop, 5/box									
	ES49	20mm x 10mm Loop, 5/box	R2010		909132	60-5180-104					
	ES51	13mm x 8mm Loop, 5/box				60-5180-101 *					
	ES52	10mm x 8mm Loop, 5/box								31-5054	
	ES53	15mm x 12mm Loop, 5/box						E1560			
	LLETZ Square Electrodes	ES14	5mm x 5mm Square, 5/box					E1565	6153901		
		ES15	10mm x 4mm Square, 5/box	S1004	9006235			E1566	6154901	DLP-SQ2	
ES16		10mm x 8mm Square, 5/box			909017		E1567	6155901	DLP-SQ1		
ES41		5mm x 10mm Square, 5/box									
ES56		10mm x 10mm Square, 5/box									

* shows 1mm difference in size

ELECTRODE SHAPES AND SIZES...

